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AFFINIS[®] heavy body / fast heavy body • surface activated

Definition

AFFINIS is a dental impression material for use in dentistry with an automatic mixing device.

Material type

Polyvinylsiloxane, addition-type, surface activated silicone elastomer.

heavy body / fast heavy body

Base: light brown

Catalyst: caramel

– ISO 4823, Type 1, high consistency

Indications

- Tray material for the 1-Step or double mix simultaneous technique
- Tray material for primary impression in 2-Step «heavy body/wash» technique (similar to putty/wash technique.)
- Tray material for the Dual-Arch bite technique (closed mouth Dual-Arch technique)
- For full and partial denture impressions
- For all border moldings, transfer copings or «pick-ups» and implant cases (once healing has taken place; and the case is ready for impression taking.)

Important

- Check intraorally if the material is completely set before removing.
- High temperatures will speed up the setting process and low temperatures will slow it down.
- For extended working time, material can be placed in the refrigerator prior to use.

Contra-indications

Provided that the product is used as directed, there are no known contra-indications.

Side effects / Interactions

Polyvinylsiloxanes have a very good biological compatibility and up to now no harmful reactions or secondary effects on patients and/or dental personnel are known. Impression materials are intended to set in the mouth of the patient. However, they should remain in the mouth not more than twice the recommended setting time. Although they have reasonably high tear strength, care should be taken that no portion of the impression material remains in the interdental spaces or in the sulcus. Undercuts should, in certain instances, be blocked out before taking the impression.

Trays

Tray selection depends on impression technique and professional preference of stock trays, Dual-Arch bite trays or custom trays. For perfect adhesion, we recommend applying a thin layer of Coltène[®] Adhesive (dry for minimum 1 min) or any other brand of adhesive specified for use with polyvinylsiloxane impression materials.

Dispensing preparation

1. Pull off the cartridge plug and discard it.
2. Eject a small amount of material directly out of the cartridge onto a paper towel/tissue held in your hand. Extrude enough material until it is evident that base and catalyst are being ejected. This ensures proper mixing.
3. Immediately wipe the cartridge orifice to remove excess material.
4. Attach the selected mixing tip to the cartridge and twist a ¼ turn clock-wise (90 °) to lock it in place.
5. Squeeze the trigger with a smooth, steady force to mix and dispense material. Avoid abrupt movements. Release the trigger to stop the flow.

Mixing tip

Disinfect after use by wiping the Mixing Tip with a disinfectant. Store in horizontal position. Always store the cartridge with the Mixing Tip attached. This seals the cartridge until it is used again and prevents contamination.

Only remove mixing tip directly before the device is to be used again, check that material is ejected evenly and attach new Mixing Tip.

Recommendations for use

«double mix technique»

While loading the tray, the dentist may begin to apply the wash-material around the preparation. Tray loading and application of wash-material around the preparation have to be completed at the same time. Keep the oral tip immersed in the material at all times to avoid trapping air. *Immediately place filled tray into the patient's mouth. Press for 2–3 s and hold in position until the material is completely set.*

«two step technique»

If a 2-Step «putty or heavy body/wash» technique is used, the primary impression, must be carefully cleaned with lukewarm water and dried before use. This will help guarantee good adhesion between the wash (low viscosity material) and tray material.

For extended working time, material can be placed in the refrigerator prior to use.

Important

Always wear gloves.

Some latex gloves and surfaces contaminated by these gloves (teeth, preparations, retraction cords, etc.) may interfere with the setting process of polyvinylsiloxanes. The material and surfaces where the impression will take place (teeth, preparations, retraction cords, etc.) should only come into contact with thoroughly washed and rinsed gloves (wash 15 s with a detergent, rinse with lukewarm tap water for another 15 s) or vinyl/nitrile gloves. Compounds containing eugenol or hemostatic substances can also impede perfect setting.

If (H₂O₂) hydrogen peroxide is used for disinfection, it is recommended to thoroughly rinse with lukewarm water afterwards in order to avoid bubble formation.

Check intraorally if the material is set before removing from the mouth.

Disinfection

The impression should be rinsed under running tap water after removal from the mouth. After rinsing, disinfection with a suitable commercial dental disinfectant solution will not affect the impression surface or dimensions. Acrylic trays should be protected against water absorption.

Optional

The following points should be observed when autoclaving impressions:

1. Only use autoclavable components (e.g. PRESIDENT Tray AC, Coltène Adhesive AC).
2. Rinse and clean impression thoroughly under flowing, lukewarm water.
3. Impressions can be autoclaved directly after taking the impression.
4. Autoclave in steam steriliser at 134°C/ 273°F (Prion program)

Important

When sterilising implant impressions, be sure clarify in advance with the manufacturer, whether or not the implant components (e.g. impression posts, etc.) are autoclavable.

Fabrication of models

The stone models can be poured after 30 min at the earliest. The impression remains dimensionally stable for a practically unlimited period of time (for at least 1 week). The surface tension will be reduced and pouring will be facilitated if the impression is briefly washed out with a detergent and rinsed thoroughly in lukewarm clear water afterwards. All industry-standard dental stone model materials (i.e. Fujirock Dental Stone, Hard Rock Dental Stone) can be used.

Electro plating

Compatible with all commercially available copper plating and silver plating baths.

Cleaning of trays

The impression can be removed with a blunt instrument. Soaking in a universal commercial solvent or petroleum ether will dissolve the Coltène[®] Adhesive. These solvents should only be used in a well-ventilated area. Trays can be cleaned and disinfected as usual.


Shelf life and storage

AFFINIS fulfils the intended purpose at least until the expiry date, which is displayed on the containers, in well sealed containers at temperatures 15-23 °C / 59-73 °F and 50 % relative air humidity. Impressions should be stored at normal room temperature, avoid exposure to heat and sun.

Caution

Federal law restricts this device to sale by or on the order of a dentist.

Marking

The expiry date and  number are shown on the package and cartridge.

Technical data ISO 4823

Measurements are based on 23 °C / 73 °F room temperature and 50 % relative humidity.

AFFINIS heavy body

Mixing time (15 ml):	0:10 min
Total working time:	0:40 min
Setting time:	2:00 min

AFFINIS fast heavy body

Mixing time (15 ml):	0:10 min
Total working time:	n.a.
Setting time:	1:20 min

Date of issue

04-2016

Only supplied to dentists and dental laboratories or upon their instructions.

For SDS see

www.coltene.com

AFFINIS[®] heavy body / fast heavy body • surface activated

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Coltène/Whaledent AG 🇨🇭