TURBO E MICRO **SMART CUBE**

UNIVERSAL ASPIRATORS

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TECHNOLOGICAL PROGRESS LEADS US TO THE SEARCH FOR INNOVATIVE SOLUTIONS.

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Cattani leads a revolution in the field of dental suction. In 2005, as the first worldwide, we have long moved beyond the concept of fixed speed suction and have radically changed the market scenario with the introduction of Turbo Smart, the first dental suction unit, with variable speed and constant head, which has resulted in a saving in energy consumption which can reach the 50%* if compared to fixed speed aspirators.

Our effort to achieve increasingly efficiency for the work of professionals has led us to the creation of Micro Smart, the suction unit designed for one or two chairs**, characterized by flexibility of installation and great energy savings. The great news about Micro Smart is the three preset speeds for saliva ejector, spray interceptor and surgical aspiration (only when installed with a single chair). In 2013 we concentrated our research on large centralized dental suction systems for clinics and hospitals: this brought to life our Maxi Smart for 15 chairs, which uses the latest technology to control head and flow rate while reducing running costs and plant room space requirements by half. Today we are introducing our Micro Smart and Turbo Smart Cube range of suction units, once again research and development have led to further electricity saving and sustainability.

* if compared to Turbo HP3, when working at the same head and flow rates.

** according to the local professional common habits; i.e. in Italy and Spain the machine is suggested for two chairs, whereas in France and Great Britain for only one. Please contact your local dealer for advices.

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
Micro Smart Cube

Cube suction units are smaller and quieter than their predecessors, but just as powerful, all thanks to our new technology. It would not be possible to have power and energy saving without the inverter, electronics, computerised programme and software.

These features are necessary to adjust head and flow rate, to achieve active self-protection and self-diagnosis, as well as energy saving, sustainability and a remote connection that allows monitor operation and intervention on the machine if needed.

In case of problems with the weather, with power or with improper use, the suction units switch to the active self-protection mode: when an event occurs that could damage the suction unit, it automatically switches operation mode to protect your unit with the emergency situation and once the problem has been corrected, it automatically switches back to programmed operation mode.

Cube range machines are always fitted with a noise-reducing box.

The use of smaller motors leads to electricity saving and sound pressure reduction. The new generations of aspirators are all equipped with historical memory. Operating temperatures, over currents and dangerous operational conditions can be shown on the display simply following the instructions in the manual.

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TURBO SMART AND MICRO SMART ARE INNOVATIVE MACHINES WITH 6 PATENTS.

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ENVIRONMENTAL SUSTAINABILITY

The new products have improved work performances, saved electricity and increased reliability. In fact, in the course of time the evolution in clinics has been as follows:

Turbo Smart Cube

- Until 2005 in a clinic with **4 chairs** it was necessary to use a suction unit weighing **22 kg** (Flux-Jet).
- After 2005 with Turbo Smart it's enough to use a suction unit weighing **10.5 kg** (Uni-Jet 75), thus sparing 11.5 kg of iron-aluminum and copper.
- In 2015, always for a clinic with 4 chairs it's enough to use a Turbo Smart Cube fitted with a suction unit Uni-Jet 40 weighing only **9 kg**.

13 kg of raw materials have been spared in total.

Micro Smart Cube

In 2015 we have replaced the Turbo-Jet 2 suction unit (**13 kg**) with Micro Smart fitted with a suction unit (Uni-Jet 25) weighing **6 kg** only.

7 kg of raw materials have been spared.

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Turbo Smart and Micro Smart are incredibly **reliable**: in the event of a functional problem, the **active self-protection** feature electronically adjusts the unit's operation, resolving the problem without damaging the machine.

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ACTIVE SELF-PROTECTION

Through inverter control and a computerized programme Micro and Turbo Smart Cube can adjust to cope with special situations.

- weather conditions, such as high temperature, the suction unit will not stop and will suffer no damage, it automatically decreases head for the time needed to restore the correct operating temperature;
- sudden wave of liquid from the surgery, the suction unit will not overflow but goes into active self-protection, the recirculation valve opens, suction slows down, the centrifugal separator absorbs the current spared by the suction unit and, once the wave has been dealt with, the suction unit will resume regular operation;
- sudden voltage change or power surge, within pre-set limits, the electronics protects the machine allowing operation to continue;
- any problems occurring during operation are highlighted on the display

WI-FI CONNECTION

Your unit can be connected wireless to the computer in the clinic. The dentist can check and modify the functioning of the aspirator at any time without moving from his practice. Modern technologies are indispensable to obtain new work performances such as: planning and control of the prevalence, constant flow and prevalence during operations, connection with the computer in the surgery and with the after-sale service. We especially recommend the wi-fi connection: it is at no extra charge and can assist in preventing the unit from suffering due to the user unaware of malfunctions and continuing to use the unit.

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
Turbo Smart Cube
versione "B" / "B" version

Turbo Smart Cube is equipped with Uni-Jet 40 suction unit replacing the Uni-Jet 75 and preserves all the features of its immediate predecessor, with a maximum flow of 1400 l/min ("A" version) and of 1700 l/min ("B" version), maximum operating head of 210 mbar, from 1 to 4 chairs in simultaneous suction.

High performance: Turbo Smart Cube can be installed in parallel without peripheral units and without connectors. It is possible to upgrade the Turbo Smart Cube model "A" (for 2 surgeries) to model "B" (for 4 surgeries) by entering a password into the machine, either before or anytime after the installation of the unit. The operator can regulate the head value according to his variable needs.

Consequently, with Turbo Smart Cube the dealer can count on one single aspirator adapting to all customers' requirements, while the dentist is sure to receive the aspirator which meets his needs, therefore there is no problem of choice for either of them.

Result: A significant advantage for everyone!

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Micro Smart Cube is equipped with a new suction unit Uni-Jet 25 and, like our Micro Smart, is recommended also for two surgeries operating simultaneously and on continuous service, with maximum operating head of 210 mbar at a maximum flow of 1000 l/min.

The flow rate of the suction unit is sufficient for 1 or 2 chairs* working simultaneously. During operation the head remains constant whether only one or both operators are working, and any minor variation shall not be noticed.

High performance: Micro Smart Cube can be installed in parallel without peripheral units and without connectors. The operator can directly regulate the head value.

* according to the local professional common habits; i.e. in Italy and Spain the machine is suggested for two chairs, whereas in France and Great Britain for only one. Please contact your local dealer for advices.

** Sound pressure level tested according to the standard ISO 3746-1979 (E). Parameters: r= 1.5 - background noise: 34 dB (A) - instrument Bruel & Kjaer type 2232.

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MICRO AND TURBO SMART CUBE ARE MORE COMPACT,
MORE SILENT, MORE EFFICIENT AND HIGH TECH.

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WE HAVE BEEN SPECIALISING WITH AIR TECHNOLOGY
FOR 50 YEARS:
SPECIALIZATION HAS GIVEN EXCELLENT RESULTS.

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HOW IS IT WE LEAD IN OUR FIELD, WHEN WE COST LESS THAN THE ALTERNATIVES? THIS IS HOW:

Constant research: this enables us to apply the latest technology to all of our products and solutions.

We enhance performance: electronic and information technology enable us to enhance the performance and reliability of our products.

We reduce costs: less maintenance and lower energy costs mean that we are always the most economical on a cost-benefit analysis.

We reduce environmental impact: we save 50% on raw materials, so that you can save between 30% and 50% on electrical consumption.