Lab-Putty

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Definition
Coltène® Lab-Putty is an elastomeric putty for use in the dental laboratory.

Material type
Polyisoxlane, condensation-type elastomer.


Indications
Coltène® Lab-Putty is exclusively designed for extraoral application. It has been conceived for the manufacture of:

- Matrix
- Repro, general occlusion and functional pattern (gfp), matrix as positioner for clasps, denture teeth (i.e. partial dentures, temporary partials, orthodontic retainers etc.), matrix for making gingival masks Coltène® Gi-Mask.
- Repair models (lower dentures with deep undercut)
- Block-outs
- Hard palate relief stamp (mould)
- Labial/buccal positioner matrix

Contra-indications
Provided that the product is used as directed, there are no known contra-indications.

Side effects / Interactions
Polyisoxlanes have a very good biological compatibility. Up to now no harmful reactions or secondary effects on laboratory staff are known.


Dosage

Measure out Coltène® Lab-Putty Base with the measuring scoop (take off excess). Spread measured out quantity on the palm of the hand. For each scoop Base press the scoop-rim once into the mass. Dose one strand-length of Coltène® Lab-Putty Activator per circle diameter.

Close tube immediately after use.

Mixing

Fold material and knead with fingertips for 30–45 s until a homogeneous, pink-colored mixture is produced.

Setting time

Min. 5 – 7 min at room temperature.
Intensive mixing, high temperatures or over-dosage of the activator speed up the setting.

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Low temperatures and under-dosage of the activator slow down the setting.
Max. over-dosage / under-dosage: 20%


Shelf life and storage
Coltène® Lab-Putty fulfills the intended purpose at least until the expiry date, which is displayed on the packaging, when stored at 15–23 °C/59–73 °F and 50% relative air humidity.

Marking
The expiry date and ☞ number are shown on the package and cartridge.

Technical data
Mixing time: 30 – 45 s
Total working time: ca. 90 s
Setting time: 5 – 7 min
Shore A hardness: 63 – 70
Mixing ration Putty : Activator 11,0 g : 0,3 g

Mixing conditions
The stated working time is based on measurements at 23 °C/73 °F room temperature and 50% relative humidity.

Lab Putty Activator – Precautions

GH507 GH508

Warning
H315 Causes skin irritation.
H319 Causes serious eye irritation.
H373 May cause damage to organs through prolonged or repeated exposure.

H413 May cause long lasting harmful effects to aquatic life.

P280 Wear protective gloves / protective clothing / eye protection / face protection.

P302+P352 IF ON SKIN: Wash with plenty of water.

P332+P313 IF skin irritation occurs: Get medical advice / attention.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313 IF eye irritation persists: Get medical advice / attention.

P314 Get medical advice / attention if you feel unwell.

P273 Avoid release to the environment.

Date of issue
06-2018

Only supplied to dentists and dental laboratories or upon their instructions.

Caution:
Federal law restricts this device to sale by or on the order of a dentist.

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