# ROEKO Gelatamp

Instructions for use
----------------------

**Proprietés:**

Gelatamp is a haemostatically acting gelatine sponge to which 5% colloidal silver is added. It facilitates optimum wound treatment when applied to a surgical cavity and can be cut to the required size to fit smaller wound cavities. The evenly porous foam structure absorbs its own weight in blood several times over, promotes thrombocyte aggregation due to the large surface and fills the wound cavity. The plug thus formed has a constant volume, fits snugly and stabilises blood coagulum. This prevents the formation of fissures and secondary cavities which, without Gelatamp, could form by contraction of the blood coagulum and trigger infection due to the invasion of contaminated saliva. Callus formation is not hindered in this way. Impaired wound healing in larger surgical cavities is thus avoided.

Gelatamp remains in the wound and is completely absorbed within four weeks.

The addition of colloidal silver has an antimicrobial effect and does not develop resistance. Unlike other potential antimicrobial additives, colloidal silver cannot be washed away from the sponge so that its insolubility produces a long lasting depot effect. Gelatamp is sold gamma-sterilised.

**Composition:**

One Gelatamp gelatine sponge (14 x 7 x 7 mm) contains:

Hardened gelatine Ph. Eur.	9,5 mg
Colloid silver Ph. Eur.	0,5 mg

**Indications:**

The treatment of alveoli and wound cavities, e.g. after extractions, cystectomies, cystostomies, apical amputations, maxillary sinus perforations, following surgical removal of tumours or retained teeth
Prevention of secondary cavity formation
Prophylaxis of wound infections
Secondary haemorrhage prophylaxis
As a dressing after gingivectomy
In periodontopathies

**Instructions for use:**

After opening the aluminium seal, Gelatamp is removed under aseptic conditions and is immediately ready for use. The size of the small sponge can be adjusted to fit the wound cavity if need be, without squeezing. Two Gelatamp plugs can be used for larger wounds.

The sponge is immediately placed dry and un-squeezed in the fresh, blood-filled wound cavity and should be allowed to absorb as much blood as possible (visual inspection). The sponge must not protrude over the inner gingival margin so as not to adversely affect the desired epithelialisation. This constitutes definitive treatment of the alveolus and there is no need for subsequent rinsing.

**Sterility:**

Gelatamp is supplied gamma-sterilised. The sterility of undamaged, unopened packs is guaranteed until the first application. Damaged packages should be discarded prior to use.

Given the antimicrobial effect of this product, safe application is guaranteed until the product is used or up to the expiry date once the package had been opened, provided that Gelatamp is removed under aseptic conditions.

**Contra-indications:**

Infected, secretory and purulent wounds; hypersensitivity to silver or gelatine.

**Side effects:**

Recovery may be delayed if used incorrectly, especially in the case of contaminated alveolar spaces.

**Interactions with other drugs:**

Concomitant application of alveolar pastes, wound plugs or other inserts apart from thrombin preparations as dry powder, and moistening with liquid substances or medication such as anion-active surfactants, antiseptics containing phenol, ethacridine lactate, dilute mineral acids and concentrated saline solutions can hinder wound recovery and should therefore be avoided. Please avoid contact with hydrogen peroxide as this will bleach the sponge.

**Storage instructions:**

Gelatamp possesses moisture-absorbing properties. The package must therefore be closed immediately after use and stored dry, protected from light. Keep all medicines out of the reach of children. Gelatamp has a shelf-life of 3 years. It can be used unopened and after opening up to the expiry date. The product should no longer be used once the date on the packaging and label has elapsed. Unused sponges can easily be disposed of with laboratory waste whereas, the plastic container, lids and aluminium foil, etc. can be recycled.

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