- (en) Automatic Mixing Unit
- (de) Automatisches Mischgerät
- (fr) Mélangeur automatique
- (it) Mescolatore automatico
- (es) Mezcladora automática de materiales de impresión
- (pt) Misturador automático
- (nl) Automatisch mengapparaat
- (el) Συσκευή αυτόματης ανάμιξης
- (sv) Automatisk blandningsapparat
- (fi) Automaattinen sekoituslaite
- (da) Automatisk blandeapparat
- (no) Automatisk blandeapparat
- (ko) 자동 믹스 유닛
- (zh) 自動混合器
- (th) เครื่องผสมวัสดุพิมพ์ปากอัตโนมัติ

Operating Instructions  
 Betriebsanleitung  
 Mode d'emploi  
 Istruzioni per l'uso  
 Instrucciones de uso  
 Manual de instruções  
 Bedieningshandleiding  
 Οδηγίες χρήσεως

Driftsanvisning  
 Käyttöohjeet  
 Betjeningsvejledning  
 Driftsveiledning  
 사용설명서  
 使用説明書  
 ขอนแนะนำในการใช้งาน

3M Deutschland GmbH  
 Health Care Business  
 Carl-Schurz-Str. 1  
 41453 Neuss – Germany


3M, ESPE, Penta and Pentamix are trademarks of 3M or 3M Deutschland GmbH.  
 Used under license in Canada. © 2020, 3M. All rights reserved.

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
3M, ESPE, Penta and Pentamix are trademarks of 3M or 3M Deutschland GmbH.  
Used under license in Canada. © 2020, 3M. All rights reserved.

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### Operating Instructions

## Table of Contents

	Page
1. Safety Instructions	5
2. Product Description	5
3. Technical Data	6
4. Installation	6
5. Operation	6
5.1 Swiveling the Unit into Working Position	6
5.2 Loading the Cartridge	6
5.3 Inserting the Cartridge	6
5.4 Attaching the Mixing Tip	6
5.5 Starting a New Cartridge	6
5.6 Loading the Impression Tray and Penta Elastomer Syringe	6
5.7 Removing the Cartridge	7
5.8 Re-starting a Partly Empty Cartridge	7
5.9 Replacing a Cartridge	7
5.10 Storage of Partially Used Cartridges	7
5.11 Fill Level Gauge	7
5.12 Parking the Unit	7
5.13 Transporting the Unit	7
5.14 Disposal of Consumable Materials	7
5.15 Return of Old Electric and Electronic Equipment for Disposal	7
5.16 Notes	7
6. Malfunctions and Error Messages	8
7. Maintenance and Care	8
7.1 Replacement of the Plunger Discs	8
7.2 Care	8
7.3 Maintenance and Repair	8
8. Customer Information	9
8.1 Warranty	9
8.2 Limitation of Liability	9
Symbol Glossary	9

## 1. Safety Instructions

ATTENTION! Read these Operating Instructions carefully before plugging in and operating the unit! As with all technical units, the trouble-free function and safe operation of this unit can be assured only if standard safety procedures and the special safety instructions in these Operating Instructions are observed.

**3M Deutschland GmbH guarantees trouble-free operation of the Pentamix system only if 3M system components are used exclusively.**

1. The unit may be used only in strict accordance with the following instructions. We accept no liability for any damage arising from any other or improper use of this unit.
2. If the unit is moved from a cold into a warm environment, condensation can result in a hazardous situation. Never begin operating the unit until it has reached the ambient temperature.
3. The unit must be placed on a non-slip, level surface to insure stable placement.
4. Always turn off the unit and disconnect the plug from the power source before changing the plunger discs.
5. When changing the plunger discs, position the smooth side of the new discs towards the foil bag (the opposite side of the discs with the fan-shaped depressions must face away). Otherwise, the foil bag might be damaged.
6. Use only genuine 3M parts when replacing defective components as directed in these Operating Instructions. We assume no liability for any damage caused by the use of non-3M parts.
7. If, for any reason, it is possible that the safety of a unit has been compromised, that unit must be removed from operation and identified in such a way that it is not inadvertently operated again by a third party. Safety may have been compromised, for example, if the unit does not work as specified or is visibly damaged.
8. Keep solvents, flammable liquids, and powerful heat sources away from the unit as they may damage the plastic housing.
9. Do not allow any liquids to enter the unit during cleaning as they could cause an electrical short or a dangerous malfunction.
10. Only authorized personnel may open the unit housing and repair the unit.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

## 2. Product Description

Pentamix™ Lite is an automatic mixing unit for the preparation of Penta™ impression materials. 3M Penta impression materials have been specifically designed for this unit and are supplied in special foil bags. This is the only form of packaging which can be used in the Pentamix Lite units. Every Penta™ foil bag is sealed with a **PentaMatic** sealing cap. The external geometry of the sealing cap insures its proper positioning on the cartridge. When the Pentamix unit is turned on, a special internal device opens the foil bag automatically once sufficient pressure has been built up by the plunger. The connection pieces located on the **PentaMatic** sealing caps serve for proper attachment of the mixing tip.

The cartridges and the **PentaMatic** sealing caps act together as a stable unit for the preparation of pastes in the Pentamix unit and allow a change of products without intermediate cleaning of the device and with no loss of paste. Do not operate the Pentamix Lite unit using any elements other than Penta™ cartridges for Pentamix Lite and Penta™ Mixing Tips – Red.

Thanks to the integrated swivel mechanism and carrying handle, the Pentamix Lite mixing unit can be folded into a compact bundle for simple transport.

These Instructions for Use should be kept for the duration of product use. The product may only be used when the product labeling is clearly readable. For details on all additionally mentioned products please refer to the corresponding Instructions for Use.

#### Intended Purpose

Intended purpose: dental impression material mixer.

Intended users: educated dental professionals, e.g. general dentists, dental assistants, dental hygienists, who have theoretical and practical knowledge on usage of dental products.

#### Contraindications

None

#### Potential Undesirable Side Effects and Complications

None known.

## 3. Technical Data

Operating mode:	Short term operation – max. 1.5 min. on, min. 10 min. off
Cartridges:	Suitable for use with Penta foil bags Dimensions: 60 x 130 mm, 26.8 x 130 mm
Housing:	All parts of the housing are made of impact-resistant plastic
Power supply:	100–120 V / 230 V (see rating plate)
Frequency:	60 Hz / 50 Hz (see rating plate)
Power input:	2.5 A / 1.5 A
Ambient temperatures:	18° C/64° F..... 40° C/104° F
Rel. humidity:	20%.....80%
Classification:	Hazard Class I
Dimensions:	L 395 mm x W 265 mm x H 305 mm
Weight:	approx. 6.7 kg
Storage and transport:	–20° C to +60° C/–4° F to 140° F Maximum rel. humidity 80%.
Emission sound pressure level:	approx. 65 dB(A).

## 4. Installation

Prior to installing the unit, make sure the power supply corresponds to the supply quoted on the rating plate of the unit.

In order to work with the unit, the power plug must be connected into a suitably grounded socket. In case of using an extension cable,  make sure that the ground line is not interrupted.

- Pentamix Lite can be placed on any non-skid flat surface. The unit is ready for operation as soon as the power cable has been plugged into an electrical outlet and the power switch has been turned on.
- The Pentamix Lite mixing unit can be permanently left in operational mode. Power consumption in the standby mode is extremely low (<1 W).

## 5. Operation

### 5.1 Swiveling the Unit into Working Position

Safe and correct operation of the Pentamix Lite is assured only when it is in working position. When swiveling the unit, use caution to avoid  the possibility of injury from jamming or crushing your fingers.

- ▶ With one hand, slowly tip the mixing apparatus of the Pentamix Lite mixing unit forward from the horizontal park position until it audibly snaps into the working position (Fig. 1).

- ▶ Do not push the unit beyond the final position as this could lead to technical problems.

### 5.2 Loading the Cartridge

Please carry out the steps in the following sequence exactly.

- ▶ Have ready one foil bag of base paste and one of catalyst paste as well as a matching 3M Penta cartridge for Pentamix Lite.
- ▶ Insert the foil bags into the unlocked cartridge (Fig. 2).
- ▶ The locating lugs make sure they are correctly seated.
- ▶ Press down on the cartridge locking lever (Fig. 3) until it audibly clicks into place.

### 5.3 Inserting the Cartridge

- ▶ Take the cartridge in your left hand.
  - ▶ Turn the plunger adjustment knob clockwise with your right hand until it stops and hold it in place. The plungers are now in the uppermost position.
  - ▶ Place the cartridge into the unit.
  - ▶ Rotate the plunger adjustment knob anticlockwise until the plungers markedly touch the foil bags (Fig. 5).
- The unit cannot be activated until the plungers touch the foil bags.
- ▶ Attach a mixing tip prior to begin mixing (see Section 5.4).

### 5.4 Attaching the Mixing Tip.

#### Always Attach a New, Unused Mixing Tip Before the Application!

Foil bags equipped with a **PentaMatic** sealing cap may be used exclusively with 3M Penta mixing tips.

- ▶ Release the lock (Fig. 6).
- ▶ First, place the mixing tip only on the tip of the mixing shaft (Fig. 7). Turn slightly to align the mixing tip with the corresponding outflow openings of the foil bags (Fig. 8). Attach the mixing tip and press down firmly to insure that the mixing tip is seated tightly on the outflow openings (Fig. 9).
- ▶ Press the locking lever down until it audibly clicks into place (Fig. 10).

### 5.5 Starting a New Cartridge

- ▶ Press the green "Start" button on the cartridge (Fig. 11).
- ▶ When a new pair of foil bags is started, you must wait **15–25 sec.** until pressure has built up in the cartridge; the foil bags will then open automatically.

Once the bags are open, one of the two pastes, which are of different colors, may flow out first, producing an uneven mixture. For this reason, **observe the flow of the two pastes, discard any uneven mixture, and do not apply the paste until the mixture shows a homogeneous color** (Fig. 12).

- ▶ After mixing, leave the mixing tip on the cartridge as a seal!

### 5.6 Loading the Impression Tray and Penta Elastomer Syringe

- ▶ Press the green "Start" button on the cartridge.  
Please note: A few seconds pass before the paste becomes visible in the mixing tip.

The unit only functions as long as the "Start" button is held down. As soon as the button is released, the paste does not flow any more. When the "Start" button is released, the Pentamix Lite unit automatically relieves the tension on the foil bags by withdrawing the plungers, which is indicated audibly by the sound of the drive motor.

- ▶ Tilt the tray to load it. Fill the tray bottom continually from one end to the other, keeping the opening of the mixing tip immersed in the impression material (Fig. 13). If necessary, smooth the impression material with a spatula.

A mixing tip may only be used as long as the mixing shaft inside the mixer rotates smoothly. If the mixing procedure is interrupted for more than 30 sec, the material inside the mixing tip begins to harden and a new mixing tip must be assembled (refer to 5.4 "Attaching the Mixing Tip"). If the mixer is operated while containing hardened material, the catalyst nipple of the sealing cap can be damaged or may break away.

- ▶ Put the Penta Elastomer Syringe (after removing the plunger) with its rear opening directly onto the mixing tip and press the "Start" button (Fig. 14).

CAUTION: Not all Penta impression materials can be syringed. Please refer to the instructions for use of the respective product.

#### ▶ **Cleaning the Penta Elastomer Syringe**

Please refer to the pertinent Instructions for Use for information about cleaning and disinfecting the Penta Elastomer Syringe.

### 5.7 Removing the Cartridge

**Do not press the "Start" button while removing the cartridge!**

- ▶ Turn the plunger adjustment knob anticlockwise with your right hand until it stops and hold it in place. The plungers are now in the upper-most position.
- ▶ Using your left hand, remove the cartridge from the device (Fig. 15).

### 5.8 Re-starting a Partly Empty Cartridge

- ▶ After prolonged storage (especially for silicones), use an instrument to remove any plugs that may have formed from the outflow openings of the foil bags.
- ▶ It may be necessary to start the unit without a mixer to check the paste flow.
- ▶ Attach a new mixer and start the mixing procedure.

### 5.9 Replacing a Cartridge

If a cartridge is nearly empty, we recommend that you have a second cartridge ready for use at hand so that you will be able to replace the exhausted cartridge quickly. This cartridge should be prepared as described under 5.5 "Starting a New Cartridge".

Do not forget to move the plungers into contact with the bags!

Lifetime of the mixing tip: see 5.6.

### 5.10 Storage of Partially Used Cartridges

- ▶ Foil bags which have been started and partially emptied must be stored only in the cartridge and with attached and filled mixer.
- ▶ Store any cartridges with partially used foil bags either horizontally or with the mixers pointing down (Fig. 16). If stored improperly, the impression material may flow backwards, which would prolong the restarting time or cause the base and catalyst pastes to mix.
- ▶ When there is a partially used cartridge in the unit and it is in park position, the plunger discs hold the foil bags in the cartridge. It is not necessary to remove the cartridge before parking the unit.

### 5.11 Fill Level Gauge

The fill level gauge shows how much material remains in the foil bags (Fig. 17).

When the foil bags are completely empty, the red fill gauge marker disappears, and the unit shuts off automatically.

### 5.12 Parking the Unit

The Pentamix Lite unit can easily be folded into a compact park position when its operation is temporarily not required.

When swiveling the unit, use caution to avoid the possibility of injury from jamming or crushing your fingers.

- ▶ Before tipping the mixing apparatus of the unit containing a cartridge, always be sure that the plunger discs are in the cartridge and that the mixing tip which has been used is still in place. The cartridge may otherwise fall out of the unit and be damaged.
- ▶ With one hand, slowly tip the mixing apparatus of the Pentamix Lite mixing unit backward from the working position until it audibly snaps into the park position (Fig. 18).
- ▶ Do not push the unit beyond the final position as this could lead to technical problems.

### 5.13 Transporting the Unit

When in the park position, the Pentamix Lite unit can be easily transported and stored.

When transporting the unit, use caution to avoid the possibility of injury from jamming or crushing your fingers.

- ▶ Move the unit into the park position before transport; see 5.12.
- ▶ Turn off the power switch on the Pentamix Lite unit and remove the power cable from the unit.
- ▶ Grasp the unit from the side with one hand (Fig. 20) and carry it using the carrying handle in the housing (Fig.19).

### 5.14 Disposal of Consumable Materials

Dispose of used mixing tips and empty foil bags in residual waste. Be sure to comply with local waste disposal regulations.

The used-off-packaging of the Penta base paste is subject to the EU Waste Electrical and Electronic Equipment Directive (WEEE) and must be collected and disposed in the EU separately from other waste in accordance with the local regulations for electric or electronic waste.

### 5.15 Return of Old Electric and Electronic Equipment for Disposal

#### 1. Collection

Users of electric and electronic equipment are required to collect their old equipment separately from other waste in accordance with the regulations of the specific country. Old electric and electronic equipment must not be disposed of with unsorted household waste. This separate collection is a prerequisite for recycling and reprocessing as an important method for preserving environmental resources.

#### 2. Return and Collection Systems

When your Pentamix Lite is no longer usable, do not dispose of the device with household waste. 3M Deutschland GmbH has set up special disposal facilities to handle the equipment. Details about the procedure for the specific country can be obtained from the pertinent 3M subsidiary.

### 5.16 Notes

Before using the foil bag, record the batch data printed on the bag (batch code, use-by date, 14-digit GTIN) for the office's own documentation or for questions to the manufacturer and archive the data so that they can later be accessed at any time.

Penta Mixing Tips - Red are single-use devices and would be damaged when removing set impression material.

Please keep the single-use accessory in the original bag to guarantee the traceability of the batch.

## 6. Malfunctions and Error Messages

Malfunction	Cause	Remedy
The unit does not start	The plungers are either fully advanced or retracted	Turn the plunger adjustment knob to move the plungers from their upper or lower limit
Unit starts, but no paste is extruded	Material has set in the mixing tip	Attach a new mixing tip
	The plungers are jammed	Turn the plunger adjustment knob to check for smooth movement
The paste flows too slowly or not at all	The pastes are too cold, or sealing cap of the base paste is defective	Allow the pastes to attain room temperature (at least 18° C). Replace defective foil bag
The paste is extruded but the rotor in the tip does not rotate	The drive shaft has not engaged the mixing tip, it is stuck in its top position	Remove the mixing tip. Turn the plungers several times up to their upper limit, thus regaining axial moveability. Attach new mixing tip
The cartridge will not fit into the unit	The plungers are not fully retracted	Turn the plungers up to their limit by turning the adjustment knob clockwise and arresting it. Check the plunger discs for damage
	The cartridge receptacle is dirty	Clean the cartridge receptacle
The unit does not switch off	The "Start" button is stuck	Press the "Start" button again or unplug the unit and release the "Start" button
The plungers are difficult to retract	Foil of a bag is wedged between the cartridge and the plunger	Retract plunger against the resistance, but avoid damage of plunger discs

## 7. Maintenance and Care

### 7.1 Replacement of the Plunger Discs

If damaged, the plunger discs must be replaced. Only original parts may be used as replacements to insure safe operation.

- ▶ Turn off the unit and disconnect the plug from the power source. Use an Allen wrench to remove the screw in the center of the plunger disc. Remove the disc and screw a new one into place (Fig. 21).

When changing the plunger discs, make sure you position them with the correct alignment (see Fig. 21) to prevent damage to the foil bags.

### 7.2 Care

Clean all components with a soft cloth and, if necessary, a mild cleaning agent (e.g., dish-washing detergent). Solvents or abrasive cleaners may not be used in any case, as they can damage the plastic components!

- Cleaning agents must not enter the unit!
- Disinfectant agents must not enter the unit!

### Cleaning and Disinfection

#### Caution

Make sure that the cleaning and disinfection agents you have chosen do not contain any of the following materials:

- Oxidizing agents (e.g., hydrogen peroxides)
- Oils
- Glutaraldehydes

Use only validated methods for cleaning and disinfection.

Use solely cleaning and disinfection agents tested for efficacy and compatibility with the medical device (e.g., FDA certification). CaviWipes™ have been tested for use as a cleaning and disinfection agent. Always observe all applicable legal and hygiene regulations for dental offices and/or hospitals.

#### General

The following procedure should be applied after use or before any subsequent use. The device must first be cleaned manually and, after cleaning, be manually disinfected in accordance with the following instructions.

#### Manual Cleaning by Wiping

- ▶ Clean the device after every application with ready-to-use disinfectant wipes (e.g., CaviWipes™ based on the active substances alcohols and quaternary compounds) or equivalent cleaning agents. Follow the cleaning agent manufacturer's directions for use.
- ▶ Thoroughly clean the material surface of the device with two ready-to-use CaviWipes™ for at least 30 seconds or until no more dirt is visible.
- ▶ Special care must be exercised in disinfecting the grooves and connecting points.

#### Manual Disinfection by Wiping (Intermediate Level Disinfection)

- ▶ Disinfect the device with ready-to-use disinfectant wipes (e.g., CaviWipes™ based on the active substances alcohols and quaternary compounds) or equivalent disinfection agents. Follow the cleaning agent manufacturer's directions for use, observing in particular the contact times.
- ▶ Disinfect the surface of the device for three minutes using at least two new disinfectant wipes. Keep all surfaces moist during the three-minute period.
- ▶ Special care must be exercised in disinfecting the grooves and connecting points.

#### Inspection, Maintenance, Testing

- ▶ Examine the device for damage, discoloration, and contamination before every use.
- ▶ Do not use damaged devices under any circumstances. If the device is not visibly clean, repeat the cleaning and disinfection procedure.

### 7.3 Maintenance and Repair

No components inside the unit require servicing or repair by the user. 3M Deutschland GmbH will accept responsibility for the safety, reliability, and correct function of the unit only if:

- Any modifications or repairs are carried out only by personnel authorized to do so by 3M Deutschland GmbH;

- The electrical fittings in the room comply with all regulatory requirements;
- The unit is operated according to these Operating Instructions.

Plastic cartridges are subject to wear depending on the material and frequency of use and are therefore not covered by the warranty for this unit. They should be regularly checked for signs of wear (e.g., fissures) and replaced with new Penta cartridges for Pentamix Lite two years after initial operation at the very latest. The use of worn or defective cartridges can have adverse effects on the evenness of the mixture, reducing the quality of the impression.

## 8. Customer Information

No person is authorized to provide any information that deviates from the information provided in this instruction sheet.

### 8.1 Warranty

3M Deutschland GmbH warrants this product to be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

### 8.2 Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

## Symbol Glossary

Reference Number and Symbol Title	Symbol	Description of Symbol
ISO 15223-1 5.1.1 Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
ISO 15223-1 5.1.3 Date of Manufacture		Indicates the date when the medical device was manufactured.
ISO 15223-1 5.1.5 Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified.
ISO 15223-1 5.1.6 Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1 5.1.7 Serial number		Indicates the manufacturer's serial number so that a specific medical device can be identified.
ISO 15223-1 5.3.1 Fragile, handle with care		Indicates a medical device that can be broken or damaged if not handled carefully.
ISO 15223-1 5.3.4 Keep dry		Indicates a medical device that needs to be protected from moisture.
ISO 15223-1 5.4.2 Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure.
ISO 15223-1 5.4.4 Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
ISO 7000-0623 This way up		Indicates to orientation to store and handle product.
CE Mark		Indicates conformity to European Union Medical Device Regulation or Directive.
Medical Device		Indicates the item is a medical device.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a dental professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1).
PAP 20		Indicates paper material is recyclable.
Low-density polyethylene		Indicates plastic component made of low-density polyethylene is recyclable.
Recycle electronic equipment		DO NOT throw this unit into a municipal trash bin when this unit has reached the end of its lifetime. Please recycle.
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

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- endringer eller reparasjoner bare utføres av personer som er autorisert av 3M Deutschland GmbH til slikt arbeid;