GC RELINE™

Prior to use, carefully read the instructions for use.

EN

GC RELINE CHAIRSIDE HARD DENTURE LINER

For use only by a dental professional in the recommended indications.

Recommended indications:


GC RELINE is an acrylic resin used for chairside hard denture linings and extension of denture borders. Its methyl methacrylate-free liquid reduces unpleasant odor and irritation of the gingival tissue.

Contraindications:

Do not use on patients with known sensitivity to methacrylate. In case of allergy, consult a physician.

Directions for use:

1. Roughen the area of denture to be relined. Clean and dry the denture thoroughly.
2. Coat the artificial teeth, labial and buccal surfaces of the denture with Vaseline or COE LUBRICANT using a cotton swab. Do not apply coating within 3mm (1/8 inch) of the peripheral border. For dentures with plastic teeth, apply a coating of COE LUBRICANT to the teeth to protect them from the GC RELINE material adhering to the tooth surface.
3. Coat the denture to be relined with BONDING AGENT using the enclosed brush. Air dry the solvent for approximately 10 seconds. Note: BONDING AGENT is volatile, be sure to close cap tightly after using.
4. Standard consistency is obtained by mixing the ratio of 1.8g powder. At the first graduation mark of the powder vial, and 1.0mL liquid at the first graduation mark of the dropper. Dispense the powder and liquid into the mixing cup and stir gently for 15 seconds.

5. Using a plastic spatula, apply the resin mix evenly onto the tissue contact surface of the denture.
6. Insert the denture into the mouth at approximately 1 minute from start of the mix. Instruct the patient to close lightly into occlusion. Muscle trim the borders as necessary.
7. For a full denture reline without undercuts: Remove the excess portion of the resin intraorally. Keep the denture in the mouth until curing is completed (approximately 5 minutes and 30 seconds from the start of the mix).
8. For a full denture reline with undercuts or a partial denture: Remove the denture when the material reaches a rubbery state (approximately 2' to 2'30" from the start of the mix). Trim excess and reseat denture in the mouth before it begins to generate heat. If denture is to remain out of the mouth, place it in warm water (approximately 40° to 50°C = 105° to 120°F) for 3 minutes until completely cured. Note: When undercuts or retentive areas are present, the denture should be removed from the mouth while the material is in a rubber-like state and before it begins to generate heat. Otherwise it may be difficult to disengage and remove the denture.
9. For denture border extensions: Mix material slightly thicker than standard powder/liquid ratio.
10. Remove the denture from the mouth at approximately 5 minutes and 30 seconds from the start of the mix. Finish and polish the denture in a usual manner.

Cautions:

1. Liquid and BONDING AGENT liquid are volatile and flammable, avoid sources of ignition.
2. Wash the brush with GC RELINE liquid or ethanol after applying the BONDING AGENT.
3. Shake the powder bottle lightly before use.
4. In case of contact with skin, wash thoroughly with water.
5. In case of contact with eyes, flush immediately with water and seek medical attention.

Storage conditions:

Keep dry. Excessive heat and freezing temperatures may impact product performance. For optimal performance store in a cool, dark place (2°C to 37°C / 36°F to 98°F). (Shelf life: 3 years from date of manufacture).

Package:

Powder
Liquid
Bonding Agent
Mixing Cups
Powder Measure
Liquid Measure

Last revised: 02/2014

  0086
MANUFACTURED BY

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803, U.S.A.

 DISTRIBUTED BY
EU: GC EUROPE N.V.
Interleuvenlaan 33, 3001 Leuven, Belgium
GC ASIA DENTAL PTE. LTD.
11 Tampines Concourse, #03-05, Singapore 528729
TEL: +65 6546 7588

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GC RELINE™

GC RELINE CHAIRSIDE HARD DENTURE LINER
For use only by a dental professional in the recommended indications.

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
Directions for use:

- Roughen the area of denture to be relined. Clean and dry the denture thoroughly.
- Coat the artificial teeth, labial and buccal surfaces of the denture with Vaseline or COE LUBRICANT using a cotton swab. Do not apply coating within 3mm (1/8 inch) of the peripheral border. For dentures with plastic teeth, apply a coating of COE LUBRICANT to the teeth to protect them from the GC RELINE material adhering to the tooth surface.
- Coat the denture to be relined with BONDING AGENT using the enclosed brush. Air dry the solvent for approximately 10 seconds. Note: BONDING AGENT is volatile, be sure to close cap tightly after using.
- Standard consistency is obtained by mixing the ratio of 1.8g powder. At the first graduation mark of the powder vial, and 1.0ml liquid at the first graduation mark of the dropper. Dispense the powder and liquid into the mixing cup and stir gently for 15 seconds.

- Using a plastic spatula, apply the resin mix evenly onto the tissue contact surface of the denture.
- Insert the denture into the mouth at approximately 1 minute from start of the mix. Instruct the patient to close lightly into occlusion. Muscle trim the borders as necessary.
- For a full denture reline without undercuts: Remove the excess portion of the resin intraorally. Keep the denture in the mouth until curing is completed (approximately 5 minutes and 30 seconds from the start of the mix).
- For a full denture reline with undercuts or a partial denture: Remove the denture when the material reaches a rubbery state (approximately 2' to 2'30" from the start of the mix). Trim excess and reset denture in the mouth before it begins to generate heat. If denture is to remain out of the mouth, place it in warm water (approximately 40° to 50°C = 105° to 120°F) for 3 minutes until completely cured.

Note: When undercuts or retentive areas are present, the denture should be removed from the mouth while the material is in a rubber-like state and before it begins to generate heat. Otherwise it may be difficult to disengage and remove the denture.

- For denture border extensions: Mix material slightly thicker than standard powder/liquid ratio.
- Return the denture from the mouth at approximately 5 minutes and 30 seconds from the start of the mix. Finish and polish the denture in a usual manner.

Cautions

- Liquid and BONDING AGENT liquid are volatile and flammable, avoid sources of ignition.
- Wash the brush with GC RELINE liquid or ethanol after applying the BONDING AGENT.
- Shake the powder bottle lightly before use.
- In case of contact with skin, wash thoroughly with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.

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Packages:

Powder
Liquid
Bonding Agent
Mixing Cup
Powder Measure
Liquid Measure

Last revised: 02/2014

'GC'

CE 0086
MANUFACTURED BY
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803, U.S.A.

DISTRIBUTED BY
EU: GC EUROPE N.V.
Interleuvenlaan 33, 3001 Leuven, Belgium
GC ASIA DENTAL PTE. LTD.
11 Tampines Concourse, #03-05, Singapore 528729
TEL: +65 6546 7588

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