3M | Elipar™ DeepCure-S

ESPE

- ⓔ **LED Curing Light**
- ⓓ **LED-Lichtgerät**
- ⓕ **Lampe à photopolymériser LED**
- ⓖ **Lampada fotopolimerizzatrice LED**
- Ⓢ **Lámpara de fotopolimerización dental de tecnología LED**
- Ⓣ **Fotoemissor LED**
- Ⓝ **LED Uithardingslamp**
- ⓔ **Συσκευή φωτοπολυμερισμού τύπου LED**
- Ⓢ **LED-lampa**
- ⓕ **LED-valokovetin**
- ⓓ **LED polymeriseringslampe**
- Ⓝ **LED-herdelampe**

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Elipar™ DeepCure-S

LED Curing Light

Table of Contents	Page
Safety	1
Product Description	3
Fields of Application	3
Technical Data	3
Charger	3
Handpiece	3
Charger and Handpiece	4
Transport and Storage Conditions	4
Installation of the Unit	4
Factory Settings	4
Initial Steps	4
Charger	4
Light Guide/Handpiece	4
Inserting the Battery	4
Battery Charging	5
Operating Status Display of the Charger	5
Power Level Display of the Handpiece	5
Operation	5
Selection of Exposure Time	5
Activating and Deactivating the Light	6
Positioning the Light Guide	6
Removing and Inserting the Light Guide from/into the Handpiece	6
Measurement of Light Intensity	6
Recommended Cure Techniques with 3M Adhesives	6
Sleep Mode	7
Acoustical Signals – Handpiece	7
Troubleshooting	8
Maintenance and Care	9
Inserting/Removing the Battery	9
Handpiece/Battery Care	9
Reprocessing the Light Guide	9
Clean Charger, Handpiece, and Glare Shield	11
Storage of the Handpiece during Extended Periods of Non-Use	12
Return of Old Electric and Electronic Equipment for Disposal	12
Collection	12
Return and Collection Systems	12
Customer Information	12
Warranty	12
Limitation of Liability	12
Symbol Glossary	13

Safety

PLEASE NOTE!

Prior to installation and start-up of the unit, please read these instructions carefully!


As with all technical devices, the proper function and safe operation of this unit depend on the user's compliance with the standard safety procedures as well as the specific safety recommendations presented in these Operating Instructions.

1. Use of the device is restricted to trained personnel in accordance with the instructions below. The manufacturer assumes no liability for any damage resulting from the use of this unit for any other purpose.
2. Prior to start-up of the unit make sure that the operating voltage stated on the rating plate is compatible with the available mains voltage. Operation of the unit at a different voltage may damage the unit.
3. Position the unit so that the power plug is accessible at all times. The power plug is used to turn the charger on and off. To disconnect the charger from the mains, remove the power plug from the electrical outlet.
4. Use only the 3M charger which is provided with the unit. The use of any other charger can result in damage to the battery.
5. CAUTION! Elipar DeepCure-S may be operated only with the supplied battery or original 3M Elipar DeepCure-S replacement batteries. The use of other batteries can impair performance or cause damage to the device.
6. CAUTION! Avoid using this device in close proximity to other equipment or with other equipment in stacked form, as this may result in improper function. If use in the manner described before is necessary, this device and the other equipment should be observed carefully in order to ensure that they work properly.
7. CAUTION. Do not stare at source. May be harmful to the eyes. Restrict exposure to the area of the oral cavity in which clinical treatment is intended. Protect patient and user from reflection and intensive scattered light by taking the appropriate precautions, e.g., glare shields, goggles, or coverings.
8. CAUTION! As is the case for all high-intensity light-curing devices, the high light intensity is accompanied by heat generation on the exposed surface. This heat can result in irreversible damage if there is longer exposure in the proximity of the pulp or soft tissue. The exposure times given in the manufacturer's instructions must be observed exactly to avoid any such damage. Uninterrupted exposure times of the same tooth surface in excess of 20 seconds and direct contact

with oral mucosa or skin must be strictly avoided. Scientists working in this field are in agreement that the irritation caused by heat generated during light curing can be minimized by taking two simple precautions:

- Polymerization with external cooling from an air flow
 - Polymerization at intermittent intervals (e.g., 2 exposures lasting 10 seconds each instead of 1 exposure lasting 20 seconds).
9. Elipar DeepCure-S may be operated only with the supplied light guide or original 3M Elipar DeepCure-S replacement light guide. The light guide has to be seen as an applied part. The use of other light guides may result in a reduction or increase in the light intensity. The product's warranty does not cover any damage resulting from the use of third-party light guides.
 10. Condensation resulting from the unit being transferred from a cold to a warm environment may be a potential risk. Hence, the unit should be turned on only after it is completely equilibrated to ambient temperature.
 11. In order to avoid electric shock do not introduce any objects into the unit with the exception of replacement parts handled in accordance with the Operating Instructions.
 12. Use only genuine 3M parts when replacing defective components as directed in these Operating Instructions. The product's warranty does not cover any damage resulting from the use of third-party replacement parts.
 13. Should you have any reason to suspect the safety of the unit to be compromised, the unit must be taken out of operation and labeled accordingly to prevent third parties from inadvertently using a possibly defective unit. Safety may be compromised, e.g., if the unit malfunctions or is noticeably damaged.
 14. Keep solvents, flammable liquids, and sources of intense heat away from the unit as they may damage the plastic housing of the charger, the seals, or the cover on the operating buttons.
 15. Do not operate the unit in the proximity of flammable mixtures.
 16. Do not allow cleaning agents to enter the unit as this may short-circuit the unit or cause potentially dangerous malfunction.
 17. Only service centers authorized by 3M Deutschland GmbH may open the unit housing and repair the device.
 18. WARNING: this device must not be modified without the manufacturer's permission.
 19. Elipar DeepCure-S must not be used in patients, or by users, with heart pacemaker implants who have been advised to be cautious with regard to their exposure to small electrical devices.
 20. Do not use Elipar DeepCure-S in patients with a history of photobiological reactions (including individuals with urticaria solaris or erythropoietic protoporphyria) or who are currently on photosensitizing medication (including 8-methoxypsoralen or dimethylchlorotetracycline).
 21. Individuals with a history of cataract surgery may be particularly sensitive to the exposure to light and should be discouraged from Elipar DeepCure-S treatment unless adequate safety measures, such as the use of protective goggles to remove blue light, are undertaken.
 22. Individuals with a history of retinal disease should seek advice from their ophthalmologist prior to operating the unit. In operating the Elipar DeepCure-S unit, this group of individuals must take extreme care and comply with any and all safety precautions (including the use of suitable light-filtering safety goggles).
 23. This unit has been developed and tested in accordance with the relevant EMC regulations and standards. It is in conformity with legal requirements. Since various factors such as power supply, wiring, and the ambient conditions at the place of operation can affect the EMC properties of the unit, the possibility that, under unfavorable conditions, there will be EMC disruptions cannot be completely excluded. If you should notice problems in the operation of this or other units, move the unit to a different location. The EMC manufacturer's declaration and the recommended separation distances between portable and mobile RF communications equipment and the Elipar DeepCure-S unit are listed in the appendix.
 24. CAUTION! Portable high-frequency communication devices including their accessories should keep a minimum distance of 30 cm to Elipar DeepCure-S. Otherwise, the performance characteristics of the device may be reduced.
 25. Prior to each use of the unit ensure that the emitted light intensity is sufficient to safely guarantee polymerization. Use the light testing area integrated into the charger for this purpose.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Product Description

Elipar DeepCure-S is a high-performance LED light source for polymerization of dental materials. It consists of a charger and a wireless handpiece powered by a rechargeable battery. The device is a medical electrical device in accordance with IEC 60601-1 and is available as a tabletop device. Wall mounting is not possible.

In comparison with conventional light-curing devices, Elipar DeepCure-S features excellent beam collimation and a uniform beam profile, directing more of the light energy to the restoration being polymerized and producing a deep, uniform, and complete curing of the restoration.

The light source is a high-performance light diode (LED). The beam emerging from the device covers the light wavelength range of 430 to 480 nm relevant, for instance, for camphor quinone products and is suitable for use with the majority of light-curing dental materials, including materials for fillings, liners, core build-ups, fissure sealings, temporary restorations, and cements for indirect restorations.

See the manufacturer's information for the exposure time required for the specific dental material.

Settable exposure times:

- 5, 10, 15, 20 sec
- Continuous mode (120 sec)
- Tack-cure mode (1 sec)

The charger is equipped with an integrated light intensity meter for the Elipar DeepCure-S.

The unit is shipped with a light guide with 10 mm diameter. It is not permissible to use the light guides of other units.

The handpiece is equipped with a sleep mode to minimize the unit's energy consumption. The handpiece switches to sleep mode once it is placed in the charger or if left unused for approx. 5 minutes outside the charger.

☞ These Instructions for Use should be kept for the duration of product use. The product may only be used when the product labeling is clearly readable. For details on all additionally mentioned products please refer to the corresponding Instructions for Use.

Intended Purpose

Intended users: educated dental professionals, e.g. general dentists, dental assistants, dental hygienists, who have theoretical and practical knowledge on usage of dental products.

Fields of Application

- Polymerization of light-curing dental materials with photo initiator for the wavelength range 430–480 nm.
 - Though the majority of light-curing dental materials are responsive in this range of wavelengths, you may wish to contact the manufacturer of the material in question.

Technical Data

Charger

Operating voltage:	100–127 V 50/60 Hz 230 V 50/60 Hz (see rating plate for factory-set voltage)
Power input:	0.08 A (230 V) 0.16 A (100–127 V)
Dimensions:	Length: 170 mm Width: 95 mm Height: 50 mm
Weight:	650 g
Classification:	Protection class II, 

Handpiece

Power supply:	lithium-ion battery, nominal voltage 3.7 V
Utilizable wavelength range:	430–480 nm
Wavelength peak:	444–452 nm

Light intensity (between 400 and 515 nm): 1470 mW/cm² -10%/+20% (independent of battery power level)

Light emission area: 60-65 mm² (optically active)

Intermittent operation: The device has been designed solely for short-term operation. Typical operating time at room temperature (23° C): 7 min, at 40° C ambient temperature: 1 min on, 15 min off (cooling-off period)

Total exposure time with new, fully charged battery: typically 120 min

Dimensions: Diameter: 28 mm
Length: 270 mm

Weight: 250 g (incl. light guide)

Charger and Handpiece

Time to charge empty battery: approx. 1.5 h

Operating temperature: 10° C up to 40° C / 59° F up to 104° F

Relative humidity: 30% up to 75%


Atmospheric pressure: 700 hPa up to 1060 hPa

Total height with handpiece inserted in the charger: 180 mm

Year of manufacture: see rating plate

Transport and Storage Conditions:


Ambient temperature: -20° C up to +40° C / -4° F up to +104° F

Relative humidity: 30% up to 75%

Atmospheric pressure: 700 hPa up to 1060 hPa

Subject to technical modification without prior notice.

Installation of the Unit

Factory Settings

The factory settings of the unit are as follows:

- 10 sec exposure time

Initial Steps

Charger

- ▶ Please ensure first that the voltage stated on the rating plate corresponds to the existing main supply voltage. The rating plate is attached to the bottom of the unit.
- ▶ Place the charger on a level surface.
- ▶ Connect the power cable of the charger to the power supply.
 - The green LED on the charger showing that the unit is ready for operation comes on. This shows that the unit is ready for operation; please refer to the section "Operating Status Display of the Charger".


Light Guide/Handpiece



- ▶ Never place the handpiece in the charger without the battery inserted in the handpiece!
- ▶ Place the glare shield on the front of the unit.
- ▶ Autoclave the light guide prior to first use.
- ▶ Then attach the light guide firmly to the handpiece.

Inserting the Battery

- ▶ Remove the protective cap from the battery and keep with the unit packaging.
- ▶ Slowly insert the battery with the thread side towards the metal housing into the handpiece until it stops. Screw the battery clockwise by hand until the seal is pressed firmly against the metal housing. The unit will not function properly if the battery is not screwed in completely!


► If the unit does not function properly, remove the battery and reinsert it as described above.

Battery Charging


- The unit contains a powerful lithium-ion rechargeable battery. The battery does not have any memory effect and can therefore be recharged at any time by placing it in the charger (see the section describing the battery power level display on the handpiece).
- The battery can also be charged in the charger independently of the handpiece.
- Place the handpiece on the charger for approx. 1.5 hours before the first use, in order to charge the new battery completely for the first time.
 - The status LED showing the operating status of the charger flashes green during the charging operation; see below “Operating Status Display of the Charger”.

Operating Status Display of the Charger

Status LED	Operating status	
	Handpiece/battery outside the charger	Handpiece/battery in the charger
Steady green light	Charger ready for operation	Charging has been completed
Flashes green	—	Battery is being charged
Steady red light	Charge contact pins are wet	Charge contact pins are wet
Flashes alternately red and green	Malfunction in the charger	Problem during charging


Power Level Display of the Handpiece

Power level LED	Operating status	
	Handpiece outside the charger	Handpiece in the charger
Steady green light	Battery is charged, handpiece ready for operation	Not possible, handpiece is in sleep mode
Steady red light	Warning for low battery capacity, remaining capacity adequate typically for 5 10-sec cycles	Not possible, handpiece is in sleep mode
Flashes red	Battery fully discharged, exposure cycle will be completed or, if in continuous mode, stopped	Problem in charging, battery is defective or cannot be charged


Operation

Selection of Exposure Time

Exposure time options: 5, 10, 15, 20 sec, continuous mode (120 sec), tack-cure mode.

- ▶ See the instructions for use for the specific dental material when selecting the exposure time.
- ▶ The indicated exposure times assume that the light guide is held at the exact position over the material being polymerized.
- ▶ If the distance between the light guide and the restoration is increased, the exposure time must be adjusted accordingly because the light intensity weakens (see graph).

Select the exposure time by pressing the TIME button.

- The selected exposure time is indicated by the 4 green LEDs.
- Each time the button is briefly pressed, the setting advances to the next (higher) value. All 4 green LEDs will be turned on for a setting of 20 sec. Pressing the button again will turn off all of the LEDs and enable the continuous mode.
- The display advances through the available settings, if the button is kept depressed.
- While exposure is ongoing, the button for selection of the exposure time is inactive.

Activating and Deactivating the Light

- ▶ Briefly press the START button; the light will turn on.
 - The LEDs first indicate the set exposure time; all 4 LEDs light up for 20 sec. Every 5 sec, as the time runs down, the LEDs will turn off one at a time; at 15 sec remaining time, 3 LEDs will still be on, at 10 sec remaining time 2 LEDs, etc.
 - The LEDs do not come on at all in continuous mode; an acoustical signal is emitted every 10 seconds.
- ▶ If desired, the light can be turned off by pressing the START button again before the exposure time is over.
- ▶ Holding down the START button activates the tack-cure mode: the unit emits a single short light pulse

which enables the defined curing of Protemp Crown temporary restorations or a light-curing cement excess (e.g., RelyX Unicem) to enable easy removal.


Positioning the Light Guide

- ▶ Rotate the light guide into the desired position for polymerization.
- ▶ To make full use of the light intensity provided, place the light guide as close to the filling as possible. Avoid directly contacting the filling material!
 - Keep the light guide clean at all times to obtain full light intensity.
 - **Damaged light guides substantially reduce the light intensity and must be replaced immediately, sharp edges may cause serious injury!**

Removing and Inserting the Light Guide from/into the Handpiece

- ▶ The light guide has a magnetic holder. Remove the light guide from the handpiece by pulling it towards the front.
- ▶ Attach the light guide to the handpiece until it is firmly in place.

Measurement of Light Intensity

The light intensity of the Elipar DeepCure-S handpiece can only be reliably determined using the Elipar DeepCure-S charger! The circular testing area is located below the scale for the light intensity on the charger. Measurements carried out using the Elipar FreeLight, FreeLight 2 or Elipar S10 charger (predecessors of Elipar DeepCure-S), or another unit will produce incorrect results due to the different light sources and the varying arrangement of the components of the units.

Caution!

- ▶ Clean the testing area with a wet cloth.
- ▶ Without applying pressure, place the light guide on the testing area such that the light-emitting opening of the rod is level with the testing area.
- ▶ Activate the lamp by pressing the START button.
 - The number of illuminated LEDs is indicative of the measured light intensity: 5 LEDs = 100%,

4 LEDs = 90%, 3 LEDs = 80%, 2 LEDs = 70%,
1 LED = 60%.

- ▶ If the light intensity is below 100% (fewer than 5 LEDs are ON), check the light guide for contamination or defects.
- ▶ Either: clean the light guide of any contamination.
Or: clean the glare shield on the handpiece; see "Care" for instructions about both procedures.
Or: replace the defective light guide with a new light guide. Or: if the measures described above do not result in any improvement, call 3M Customer Service or the appropriate dealer.

Recommended Cure Techniques with 3M Adhesives

- Rotate the light guide into the desired position for polymerization.
- A barrier sleeve may be used. Note that sleeve can impact power and irradiance of a curing light. It is critical to ensure the sleeve is powder free and the seam is not over the tip of the light guide and that the sleeve fits well. Intermediate level disinfection (manual disinfection with ready-to-use disinfectant wipes) is indicated when a sleeve is used, otherwise the light guide should undergo high level disinfection (manual disinfection using disinfectant bath, see "Reprocessing the Light Guide").

The following sleeves have been tested:

Sleeve	Impacts Power or Irradiance
Cure Sleeve® Tip Covers (short) – Kerr	No
Complete Curing Light Sleeves-Demi – Kerr	No
TIDIShield™ Custom Fit Curing Light Sleeves – TIDI Products	No
SmartLite® Max Lens Cover Sleeve – Dentsply	No
Curelastic™ Cure-Light Sleeve (Steri-Shield)	Yes

Sleep Mode

Once the handpiece is placed in the charger, all internal functions and LEDs are automatically turned off as the handpiece switches to sleep mode. This reduces the power consumption of the battery to a minimal level. Outside the charger, the handpiece is also switched to sleep mode, if it is not used for approx. 5 min.

- ▶ To terminate the sleep mode press the START button.
 - The sleep mode termination signal (two short audible signals) is emitted indicating that the handpiece is ready for operation; the handpiece displays the latest selected exposure mode and time settings.

Acoustical Signals – Handpiece

An acoustical signal is emitted

- every time a button is pressed,
- every time the light is turned ON,
- 1 time after 5 sec exposure time, 2 times after 10 sec, 3 times after 15 sec.
Exception: in continuous mode; an acoustical signal is emitted every 10 seconds.

Two acoustical signals are emitted

- every time the sleep mode is terminated by pressing the START button,
- every time the light is turned OFF.

A 2-sec error signal is emitted, if

- the handpiece over-heats,
- the battery lacks sufficient charge.

The acoustical signals from the handpiece can be turned off (except for the 2-sec error signal). Follow these instructions to turn them off. Put the handpiece in sleep mode, e. g., by setting it in the charger. Take the device from the charger; press first the TIME button, then the START button. This changes the sleep mode of the handpiece back to active and the status "Acoustical signals activated" to "Acoustical signals deactivated". The acoustical signals can be reactivated by following the same procedure.

Troubleshooting

Error	Cause ► Solution
The charge level display on the handpiece glows red steadily.	The remaining battery charge is adequate only for typical 5 10-sec exposure cycles. ► Place the handpiece in the charger and re-charge the battery.
The power level display of the handpiece flashes red. The ongoing exposure is interrupted (light off signal is emitted) followed by a 2-sec error signal; the handpiece switches to "sleep" mode and resists further activation.	The battery lacks sufficient charge. ► Place the handpiece in the charger and re-charge the battery.
The power level display of the handpiece flashes red while the handpiece is in the charger.	Charging problem. The battery is defective or at the end of its useful life. ► Replace the battery.
The handpiece has not been used for a long time and now it cannot be turned on.	There is not enough charge in the battery to turn on the handpiece. ► Place the handpiece in the charger and re-charge the battery.
The light emission does not start when the START button is pressed; an error signal sounds for 2 sec.	The information about intermittent operation under the header "Technical Data – Handpiece" has not been observed. The handpiece has become overheated in the course of use. The handpiece can be used again once it has cooled down. ► Allow the handpiece to cool for 3 minutes, then start the next exposure by pressing the START button.

Error	Cause ► Solution
During light emission in continuous mode, an error signal sounds for 2 sec, the emission is stopped, and the handpiece changes to sleep mode.	The information about intermittent operation under the header "Technical Data – Handpiece" has not been observed. The handpiece has become overheated in the course of use. The handpiece can be used again once it has cooled down. ► Allow the handpiece to cool for 3 minutes, then start the next exposure by pressing the START button.
The LED showing the operating status in the charger glows red steadily.	The charge contact pins in the handpiece or the charger are wet. ► Dry the charge contact pins. Be careful not to bend the flexible charge contact pins in the charger.
The LED showing the operating status in the charger flashes alternately red and green. The handpiece is not in the charger.	Charger is defective. ► Have the charger repaired.
The LED showing the operating status in the charger flashes alternately red and green. The handpiece is in the charger.	Malfunction during charging. ► The battery is damaged. ► Replace the battery.
The LED showing the operating status in the charger does not come on although the plug has been connected to an electrical outlet.	Power outlet carries no voltage. ► Use a different power outlet.
	Charger is defective. ► Have the charger repaired.
The light intensity is too low.	► Clean the light guide and the protecting glass in the light guide mounting hole (please refer to "Cleaning the Light Guide").

Error	Cause ► Solution
The dental material does not cure completely.	<ul style="list-style-type: none"> ► Clean the light guide and the protecting glass in the light guide mounting hole (please refer to "Cleaning the Light Guide"). ► Check to see that the correct light guide has been attached.
The light guide cannot be attached to the handpiece.	<ul style="list-style-type: none"> ► The light guide is not designed for use with the Elipar DeepCure-S.

Maintenance and Care

The Elipar DeepCure-S unit is maintenance-free. No periodic maintenance is required. See the information contained in this chapter to secure problem-free operation.

Inserting/Removing the Battery

Never place the handpiece in the charger without the battery inserted in the handpiece!

Use 3M batteries only! The use of other manufacturer batteries or non-rechargeable batteries/primary batteries is a potential hazard and may damage the unit.

- Remove the battery from the handpiece by turning it counter-clockwise as shown by the direction of the arrow on the underside of the handpiece.
- Insert the (new) battery into the handpiece and turn it clockwise until the seal is firmly pressed against the metal housing.
- Place the handpiece in the charger for 1.5 hours in order to fully charge the new battery for the first time.
 - The LED showing the operating status on the charger flashes green; see also "Operating Status Display of the Charger".

Handpiece/Battery Care

- Use only the 3M charger which is provided with the unit; failure to do so could lead to damage to the battery!
- Do not immerse the battery in water or incinerate! Please also observe the chapter on "Safety".

Reprocessing the Light Guide

The light guide is not sterile when delivered and must be reprocessed before being used for the first time and after each patient use.

Follow manufacturer's instructions for preparation, temperatures, contact time and rinsing off the cleaning¹ and disinfecting solution² exactly.

An enzymatic detergent should be used as the cleaning solution.

Material Resistance

Make sure that the cleaning and disinfectant agents you have chosen do not contain any of the following materials:

- Organic, mineral, and oxidizing acids (minimum acceptable pH value 5.5)
- The pH value of the cleaning and disinfecting solution should be between pH 7 and 11
- Oxidation agents (e.g., hydrogen peroxide)
- Halogens (chlorine, iodine, bromide)
- Aromatic/halogenized hydrocarbons

Observe the manufacturer's information about the cleaning and disinfecting agents.

The light guide must not be exposed to temperatures higher than 135 °C (275 °F).

The light guide has been tested for up to 500 sterilization cycles.

Limitations on Reprocessing

This medical device is safe to use if it is reprocessed according to these instructions and no visible damage is seen.

Inspect the device carefully before each use and discard if there is any damage.

Point of use / Pre-Treatment

The contaminated light guide must be placed in an approved, sealed container during transport from the point of use to the reprocessing area.

The pre-treatment must be carried out before either automatic or manual cleaning and disinfecting. Immediately after using (within a maximum of 2 hours), remove gross contaminations from the light guide.

For this, rinse the light guide thoroughly under running water (drinking water quality as a minimum) for at least 10 seconds.

Use a soft brush or a soft cloth to manually remove any remaining contamination. Remove polymerized composite with alcohol; a plastic spatula may help in removing the material. Do not use any sharp or pointed tools to protect the surface of the light guide from scratching.

Cleaning and Thermal Disinfecting (Automated Washer-Disinfectors / CDD (Cleaning and Disinfection Device))

- After the light guide has been pre-treated, place it in the wire or sieve basket of the washer-disinfectant. Make sure that the devices do not touch one another in the washer-disinfectant.

- ▶ Remove the light guide from the washer-disinfector unit at the end of the drying cycle and store it under low germ count conditions.
- ▶ Examine the light guide for visible impurities. If any impurities are visible, repeat the cycle in the washer-disinfector.

Automated washer-disinfectors that meet the following requirements are suitable for use:

- The device generally has a certified effectiveness in accordance with ISO 15883. The washer-disinfector is regularly serviced and tested.
- If available, use a tested cycle for thermal disinfection (AO value > 600 or, for older devices, at least 5 min at 90 °C/194 °F).
- The water used for rinsing is minimum drinking water quality. Use low-germ and low-endotoxin water for final rinsing.
- The air used for drying should be filtered (oil-free, with low microbiological and particulate contamination).

Note: The evidence of suitability for automatic cleaning and disinfection took place using the Miele G7336 CD washer-disinfector (Gütersloh) and the cleaning agent Thermosept RKN-zym (Schülke & Mayr). The program D-V-THERMO0905, without drying time was used under worst case conditions.

Cleaning steps	Parameters	
Pre-rinse	Temperature:	10 ° ± 2 °C
	Time:	1 min
Cleaning	Cleaning temperature:	45 ° ± 2 °C
	Time:	5 min
	Cleaning solution:	Thermosept RKN-zym
	Concentration:	0.3% (3 ml/l)
Final rinse	Temperature:	10 ° ± 2 °C
	Time:	2 min
Thermal disinfection	Temperature	90 ° ± 2 °C
	Time:	5 min

Cleaning and disinfection of the light guide by wiping is only permitted in combination with a corresponding protective sleeve

Manual Cleaning

- ▶ Clean the light guide after every application for 30 seconds with ready-to-use disinfectant wipes (e.g.,

CaviWipes™ based on the active substances alcohols and quaternary compounds (QAV) or equivalent cleaning agents.

- ▶ Follow the cleaning agent manufacturer's directions for use.
- ▶ Clean with the wipes until no contamination is visible.
- ▶ Then rinse with water (drinking water quality as a minimum) and dry with a soft, lint-free cloth.

Manual Disinfection (Intermediate Level Disinfection)

- ▶ An intermediate level disinfection is possible only if the light guide is covered with a protective sleeve during use.
- ▶ Disinfect the light guide with ready-to-use disinfectant wipes (e.g., CaviWipes™ based on the active substances alcohols and quaternary compounds (QAV) or equivalent disinfection agents. Follow the disinfection agent manufacturer's directions for use, observing in particular the contact times.
- ▶ Then rinse with water (drinking water quality as a minimum) and dry with a soft, lint-free cloth.

Manual Cleaning and Disinfection the Light Guide by Immersion Bath

Manual Cleaning

- ▶ A neutral enzymatic cleaning agent is recommended (e.g., Cidezime/Enzol from Johnson & Johnson).
- ▶ Place the light guide for the specified application time into the solution, making sure that it is completely covered.
- ▶ Cleaning takes place in the immersion bath with a soft brush for at least 5 minutes.
- ▶ Then remove the light guide from the solution, thoroughly rinse with low-germ-count water (drinking water quality as a minimum), and dry with a lint-free cloth.

Manual Disinfection (High Level Disinfection)

- ▶ A disinfecting agent containing o-Phthalaldehyde (e.g., Cidex OPA by Johnson & Johnson) is recommended. It is compatible with the light-curing device.
- ▶ To disinfect, place the cleaned light guide for the specified application time (≥ 12 min) into the solution, making sure that it is completely covered.
- ▶ The instruments in the disinfection bath must not touch each other during this process.
- ▶ Then remove the light guide from the solution, thoroughly rinse with low-germ-count water (drinking water quality as a minimum) for at least 3 minutes, and dry with a lint-free cloth.
- ▶ Then check the light guide (see "Check").

¹ Cleaning agent Thermosept RKN-zym (Schülke & Mayr) for the washer-disinfector and Cidezime® /

Enzol[®] and CaviWipes[™] were used for manual cleaning.

² Cidex OPA, ready-for-use solution, was used for validation of the high level disinfection and CaviWipes[™] for validation of intermediate level disinfection.

Packaging for Sterilization

Package the light guide in single wrapped sterilization pouches (one-way packaging or double packaging).

Sterilization

Effective cleaning and disinfection are absolutely essential requirements for effective sterilization.

A fractionated vacuum procedure (at least 3 vacuum stages) or a gravity procedure (with adequate product drying) must be used for the sterilization. The use of the less effective gravity procedure is permissible solely if the fractionated vacuum procedure is unavailable. The user is solely responsible for selection of the parameters required for the gravity procedure; the user must therefore calculate or validate these parameters in accordance with the loading configuration. The drying time actually required is also directly dependent on parameters that are in the sole responsibility of the user (loading configuration and density, sterilizer condition, etc.) and must consequently be determined by the user.

- Steam sterilizer in conformity with DIN EN 13060/DIN EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- Maximum sterilization temperature 135 °C (275 °F)
- The following cycles were validated for the light guide:

	Fractional prevacuum		Gravity displacement
Temperature	132 °C (270 °F)	134 °C (273 °F)	121 °C (250 °F)
Sterilization time	4 min	3 min	20 min
Drying time	30 min	30 min	30 min

Check

Before using the light guide again, check it for damaged surfaces, discoloration, and contamination; do not use damaged light guides. If the light guide is still contaminated, repeat all reprocessing steps.

Clean and Disinfect Handpiece and Glare Shield

- Cleaning and disinfecting solutions must not enter the unit.

- Wipe the surfaces with ready-to-use disinfectant wipes³ (e.g., CaviWipes[™] based on the active substances alcohols and quaternary compounds (QAV) or equivalent disinfection agents). Subsequently inspect the surfaces for visible impurities. If any impurities are visible, repeat the cleaning procedure.
- Solvents or abrasive cleaners may not be used as they can damage the components.
- Wipe the device surfaces with a tuberculocidal disinfectant wipe³ for disinfection and keep it wet according to the disinfection agent manufacturer's directions for use, particularly observing the contact time.
- Disinfectants must not be sprayed directly onto the device surface.
- Wipe the device surfaces with a soft cloth moistened with water (drinking water quality as a minimum) and dry with a soft, lint-free cloth. If you do not rinse well with water, the components will be damaged.

Additional notes: Ensure that the charge contact pins stay dry and are not touched by any metallic or greasy parts. Do not bend the charge contact pins while they dry. Wet charge contact pins cause an operating error (malfunction message: the LED showing the operating status of the charger flashes red).

³ CaviWipes[™] were used for validation of the cleaning and disinfection procedure.

Clean Charger, Handpiece, and Glare Shield

Clean all components with a soft cloth and, if necessary, a mild cleaning agent (e.g., dish-washing detergent). Solvents or abrasive cleaners may not be used in any case, as they can damage the plastic components!

- Cleaning agents must not enter the unit!

To disinfect all components, spray the disinfectant on a towel and use it to disinfect the unit. Do not spray the disinfectant directly on the device.

- Disinfectant agents must not enter the unit!
- Dry residual disinfectants on the device with a soft and fluff-free cloth, as they damage the plastic components.

If necessary, ask the manufacturer of the disinfectant if its constant use will damage plastic surfaces.

Clean the protection glass with a soft and fluff-free cloth. Beware of scratches!

- Make sure that charge contact pins remain dry and are not contacted by metallic or greasy parts. Do not bend the charge contact pins during drying. Wet charge contact pins will cause an operating error (malfunction message: the LED showing the operating status in the charger glows red).

Storage of the Handpiece during Extended Periods of Non-Use

- ▶ If the handpiece is not to be used for an extended period of time – e.g., during vacation – fully charge the battery prior to departure or keep the handpiece inserted in the operational charger. A safety switch within the battery prevents a total discharge.
- ▶ Discharged or nearly discharged batteries must be recharged as soon as possible.

Return of Old Electric and Electronic Equipment for Disposal

Collection

Users of electric and electronic equipment are required to collect their old equipment separately from other waste in accordance with the regulations of the specific country. Old electric and electronic equipment must not be disposed of with unsorted household waste. This separate collection is a prerequisite for recycling and reprocessing as an important method for preserving environmental resources.

Return and Collection Systems

When your Elipar DeepCure-S is no longer usable, do not dispose of the device with household waste. 3M Deutschland GmbH has set up special disposal facilities to handle the equipment. Details about the procedure for the specific country can be obtained from the pertinent 3M subsidiary.

Customer Information

No person is authorized to provide any information that deviates from the information provided in this instruction sheet.



Warranty



3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

Symbol Glossary

Reference Number and Symbol Title	Symbol	Description of Symbol
ISO 15223-1 5.1.1 Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
ISO 15223-1 5.1.3 Date of Manufacture		Indicates the date when the medical device was manufactured.
ISO 15223-1 5.1.5 Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified.
ISO 15223-1 5.1.6 Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1 5.1.7 Serial number		Indicates the manufacturer's serial number so that a specific medical device can be identified.
ISO 15223-1 5.3.7 Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed.
ISO 15223-1 5.3.8 Humidity limitation		Indicates the range of humidity to which the medical device can be safely exposed.
ISO 15223-1 5.3.9 Atmospheric pressure limitation		Indicates the range of atmospheric pressure to which the medical device can be safely exposed.
ISO 15223-1 5.4.4 Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
CE Mark		Indicates conformity to European Union Medical Device Regulation or Directive.
Medical Device		Indicates the item is a medical device.
Follow instructions for use		To signify that the instructions for use must be followed.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a dental professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1).

Reference Number and Symbol Title	Symbol	Description of Symbol
Type B applied part		To identify a type B applied part complying with IEC 60601-1.
Class II equipment		To identify equipment meeting the safety requirements for Class II equipment per IEC 61140.
Use Indoors		Indicates medical device be used indoors
Regulatory Compliance Mark		Indicates product complies with applicable Australian regulatory requirements.
Battery		Battery power level
PAP 20/21		Indicates paper material is recyclable.
Recycle electronic equipment		DO NOT throw this unit into a municipal trash bin when this unit has reached the end of its lifetime. Please recycle.
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

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Informasjonsstatus: januar 2020

Appendix

Guidance and manufacturer's declaration – electromagnetic emission		
<p>The Elipar DeepCure-S Dental Light Unit is intended for use in the electromagnetic environment specified below. The customer or the user of the Elipar DeepCure-S Dental Light Unit should assure that it is used in such an environment.</p>		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Elipar DeepCure-S Dental Light Unit uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Elipar DeepCure-S Dental Light Unit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Immunity test	Compliance to IEC 60601 test levels	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Radiated RF EM fields IEC 61000-4-3	3 V/m 80 MHz - 2,7 GHz 80 % AM at 1 kHz	
Proximity fields from RF wireless communications equipment IEC 61000-4-3	see following table: Test specifications for enclosure port immunity to RF wireless communications equipment	
Electrical fast transients / bursts IEC 61000-4-4	± 2 kV for power supply lines 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Conducted disturbances induced by RF fields IEC 61000-4-6	$3 V_{rms}$ 0,15 MHz - 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	
Surges Line-to-line IEC 61000-4-5	$\pm 0,5$ kV, ± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, and short interruptions on power supply input lines IEC 61000-4-11	0 % U_T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0 % U_T ; 1 cycle and 70 % U_T ; 25130 cycles Single phase: at 0° 0 % U_T ; 2501300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Elipar DeepCure-S Dental Light Unit requires continued operation during power mains interruptions, it is recommended that the Elipar DeepCure-S Dental Light Unit be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a. c. mains voltage prior to application of the test level..		

Test specifications for enclosure port immunity to RF wireless communications equipment

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1,8	0,3	27
450	430 – 470	GMRS 460, FRS 460	FM \pm 5 kHz deviation 1 kHz sine	2	0,3	28
710 745 780	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0,3	9
810 870 930	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0,3	28
1 720 1 845 1 970	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0,3	28
2 450	2 400 – 2 570	Bluetooth, WLAN, 802 .11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0,3	28
5 240 5 500 5 785	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0,3	9

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Health Care Business

Carl-Schurz-Str. 1

41453 Neuss – Germany

3M Customer Care/MSDS Information: U.S.A. 1-800-634-2249.

3M, ESPE, Elipar, FreeLight, Protemp and RelyX are trademarks of 3M or 3M Deutschland GmbH.

Used under license in Canada. © 2020, 3M. All rights reserved.

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3M, ESPE, Elipar, FreeLight, Protemp et RelyX sont des marques de commerce de 3M ou 3M Deutschland GmbH.