

Palapress® / Palapress® vario – The classic self-curing denture acrylic.



Palapress and Palapress vario are denture acrylics for repair work, rebasing, and adding saddles to partial denture bases. Thanks to their particular handling characteristics, all tasks are completed safely, simply, and very efficiently.

The best thing: You can distinctly improve your fabrication process by utilising the different polymerisation sequences of Palapress and Palapress vario.

Palapress vario for adding up to two saddles in one session, indirect rebasing, restorations, and enlargements.



Compact Knowledge: Consistent denture colours

A “beautiful” and natural-looking denture colour is the precondition of aesthetic dentures. Pala offers a wide range of veined or unveined colours. Every batch is colour-metrically measured during production to reduce variations within the complicated mixture of pigments. This is how we guarantee batch to batch and product to product colour consistency. You achieve the best possible aesthetic results even with a combination of different materials.

It goes without saying that all Pala denture materials are cadmium free.

