

hygiene

Supra- and sub-gingival scaling



Slim tip

Interproximal spaces scaling. Finer and longer than tip No.1, it is also powerful and robust.

Supra-gingival scaling and interproximal spaces



Interproximal spaces

Its anatomical shape allow fast and efficient procedure.

periofine

Smooth biofilm elimination



Dental plaque and sub-gingival small deposits removal

Oriented tangentially: its shape adapts to the anatomy of the tooth for a painless and easy access.



Interproximal scaling of narrow areas

Left-oriented for an easy access to premolars and molars.



Interproximal scaling of narrow areas

Right-oriented for debridement and cleaning of medium pockets.

