



## elements™ Connect & Apex Connect

Cordless Endodontic Motor and Electronic Apex Locator

### **Seamless Connection. Greater Confidence.**

Add confidence to your shaping, with intuitive elements Connect and Apex Connect. Designed to work seamlessly together, interconnected system supports patient safety while simplifying your procedures.

# elements™ Connect and Apex Connect

## An Interconnected System



### Dual Light Indicators

With intuitive light indicators in both motor handpiece and the Apex file clip, you'll have greater visibility for more confident shaping.



### Patient Safety

The intuitive, easy-to-see smart light and audible indicators, as well as a large, bright viewing display alert you when you're reaching the apex. May be used independently or as a connected system.



### Portable Convenience

Easily move between operatories with this interconnected system; lightweight, ergonomic and cordless design.



## Versatile Motions

With Adaptive, Reciprocation and Rotary motions in one motor, elements Connect allows offices to maintain their existing techniques and protocols.

## Preset

9 Kerr preset file settings, 1 reciprocating setting and 4 customizable settings.

## Continuous Use



Two removable batteries are included to ensure one is always charged so you can complete your procedures safely without interruption.




## Built-in Accuracy

Ensure measurement accuracy with Apex Connect's automatic self-calibration feature that activates every time the device is turned on.

## Shape with Confidence

View real-time working length information during procedures with a smart light indicator that changes color, as well as audible indicators, as you reach the apex.



# Ordering Information

## elements Connect

815-1801	elements Connect Cordless Motor Kit
815-1906	elements Connect AC/DC Adapter
815-1908	elements Connect Contra Angle
815-1909	elements Connect Charging Base
815-1910	elements Connect Battery Pack

## Apex Connect

815-1803	Apex Connect Kit
815-1907	Apex Connect AC/DC Adapter
815-1912	Apex Connect Probe Cord
815-1913	Apex Connect Lip Hook
815-1914	Apex Connect File Holder Adapter
815-1915	Apex Connect File Holder A
815-1916	Apex Connect File Holder B
815-1917	Apex Connect Replacement Battery

## Power Cord

815-1805	elements/Apex Connect Power Cord
----------	----------------------------------



For more information scan the QR Code  
or visit [to.kerrdental.com/en/elements-connect-apex-connect](https://to.kerrdental.com/en/elements-connect-apex-connect)

**To Learn More About  
Kerr Endodontic Solutions  
Contact Your Kerr Rep**

