



Individo® Lux · Profibase

PROVEN TRAY BLANKS AND BASE PLATES FOR
PARTIAL, FULL AND IMPLANT PROSTHETICS

Individo® Lux · Profibase

A SOLID BASE IN THE LABORATORY

Individo® Lux

Indications

Fabrication of individual trays for use in crown, bridge and total prosthetic techniques

Fabrication of bases for bite patterns

Fabrication of odontograms

Advantages

- In blue-transparent or blue-opaque colour
- Allows precise and economic processing
- Stable trays
- Pleasant mint flavour
- Non-sticky consistency
- Light-cured with standard halogen light or UVA lights



Presentation

REF 2402	Maxillae 50 pcs. blue opaque
REF 2403	Mandibles 50 pcs. blue opaque
REF 2408	Maxillae 50 pcs. blue transparent
REF 2409	Mandibles 50 pcs. blue transparent

Profibase

Indications

Fabrication of individual trays for use in crown, bridge and total prosthetic techniques

Fabrication of bases for bite patterns

Base plates for wax contouring in the total prosthetic technique

Advantages

- Gingiva shade
- Simple and time-saving shaping by hand
- Permanent stability and compressive strength
- Pleasant mint flavour
- Non-sticky consistency
- Light-cured with standard halogen light or UVA lights



Presentation

REF 2111	Maxillae / mandibles 2 × 25 pcs.
REF 2112	Maxillae 50 pcs.
REF 2113	Mandibles 50 pcs.

VOCO GmbH
Anton-Flettner-Straße 1-3
27472 Cuxhaven
Germany

Tel.: +49 (0) 4721-719-0
Fax: +49 (0) 4721-719-140

info@voco.com
www.voco.com

Available from: