

ESTECEM II PLUS

The adhesive allrounder – attaches indirect restorations regardless of the material.



Colour stability

No solubility after light curing – colour stability and no wash-out of the margin.

Viscosity

Optimum viscosity – no uncontrolled flow of the paste.

Handling

Good processing – easy removal of the excess.

Shades - 4:



Find out more:



| Colour stability | | ESTECEM II | RelyX Ultimate | Multilink Automix | Variolink Esthetic DC | Panavia V5 | Panavia F2.0 |
|---|-----------|------------|----------------|-------------------|-----------------------|------------|--------------|
| Colour | Universal | A1 | Yellow | Bright | Universal | Bright | |
| Before Test | | | | | | | |
| After Test (24hours at 80 °C in coffee solution) | | | | | | | |
| ΔE | 4,6 | 8,8 | 11,0 | 19,8 | 12,5 | 11,8 | |

Source: TOKUYAMA DENTAL R&D



Customers also bought this additional product:

UNIVERSAL BOND II

More at BONDING >> UNIVERSAL BOND II